

CHEM202

Organic Chemistry

Nick Green

nick.green@otago.ac.nz, 3n11

Module 1: Spectroscopy and Characterisation


Lecture 8: NMR Spectroscopy

Certain NMR experiments can be used to
distinguish enantiomers!




Lecture 8

- Chirality and ^1H NMR
- COSY, 2D-NMR.
- NOE, NOESY

Chirality and ^1H NMR



L-aspartic acid



The CH_2 group of aspartic acid displays resonances characteristic of the AB part of an ABX system. The protons are inequivalent and are said to be **diastereotopic**.

There is free rotation about the bond, but **H1 and H2 will never occupy an equivalent pair of environments.**


The test for diastereotopic protons is: if you replace either with a deuterium atom, will it create a pair of diastereomers?


Correlation Spectroscopy (COSY)

- COSY provides a map of all coupling networks between protons in a molecule spectroscopy in a single experiment (cf with spin decoupling).
- **Coupling provides connectivity through the bonds**
- It is an example of a two-dimensional (2D) NMR technique.
- Refer Appendix 12; Advanced NMR techniques **(Expt. 3, 4 & 5)**

1D spectrum for 2-butanone

Coupling for this simple molecule can be extracted from the 1D spectrum. But what does a 2D spectrum look like?






1D spectrum on
x and y axis

‘Diagonal’ is the
same peak on
each axis


‘Cross-peaks’
appear off the
diagonal, and
only between
coupling pairs

No cross-peak
when no
coupling




COSY spectrum of 1-nitropropane.

Notes: "A" shows that "a" is coupled to "b" and "B" shows that "b" is coupled to "c."




COSY spectrum of 2-methyl-3-pentanone



A spin system is a group of protons connected (indirectly) by coupling.



This molecule has two discrete spin systems.




Iodobutane




All four signals are part of one spin system.

Thus, we expect all signals to be (indirectly) connected through the COSY

Iodobutane



**Predict the COSY
for this molecule**




All four signals are part of one spin system.

Thus, we expect all signals to be (indirectly) connected through the COSY


NOESY (Nuclear Overhauser Effect Spectroscopy)

- Enables the molecular geometry to be defined
- Effect (NOE) is transmitted **through space** (not bonds) between interacting nuclei
- Interaction is based on relaxation processes following excitation of the nuclei
- The NOE rapidly diminishes with distance (depends on d^{-6})
- Nuclei must be relatively close together ($d < 0.5$ nm)
- Refer Appendix 12; Advanced NMR techniques (**Expt. 4**)

A synthetic unknown, $C_{13}H_{10}N_2$
- could be structure A or B?




A




B

A synthetic unknown, $C_{13}H_{10}N_2$
- could be structure A or B?




A




B

**Spin systems
colourised.
The difference in
number of
protons in each
spin system will
be diagnostic.**


COSY spectrum of unknown




A synthetic unknown, $C_{13}H_{10}N_2$
- could be structure A or B?




A



B





COSY spectrum of unknown:




NOESY spectrum of unknown:

COSY spectrum of unknown





COSY spectrum of unknown:



NOESY spectrum of unknown:



Use the number of signals expected in each set of spin systems to decide.

(Note: because the COSY is always symmetrical about the diagonal, we only need one axis shown)

